Instructions for use MEYCO diamond knives (applies to ME100-ME799)

Content:

| | |
|---|----|
| (DE) Gebrauchsanweisung MEYCO Diamantskalpelle (gilt für ME100-ME799)..... | 5 |
| (DK) Brugsanvisning MEYCO diamantskalpeller (gælder for ME100-ME799)..... | 17 |
| (GB / US) Instructions for use MEYCO diamond knives (applies to ME100-ME799)..... | 2 |
| (ES) Manual de instrucciones de los bisturís de diamante MEYCO (válido para ME100-ME799)..... | 8 |
| (FR) Mode d'emploi Bistouris en diamant MEYCO (pour ME100-ME799)..... | 11 |
| (IT) Bisturi diamantati MEYCO - Istruzioni per l'uso (valide per ME100-ME799) | 14 |
| (NO) Bruksanvisning for MEYCO Diamantskalpeller (gjelder for ME100-ME799) | 20 |
| (NL) Gebruiksaanwijzing MEYCO diamantscalpels (geldt voor ME100-ME799) | 26 |
| (PL) Instrukcja użytkowania skalpeli diamentowych MEYCO (dotyczy modelu ME100-ME799) | 29 |
| (PT) Instruções de utilização dos bisturís de diamante MEYCO (válidas para ME100-ME799) | 32 |
| (SE) Bruksanvisning MEYCO diamantskalpeller (gäller för ME100-ME799)..... | 23 |

(GB / US) Instructions for use MEYCO diamond knives (applies to ME100-ME799)

Legend:

| Manufacturer | Date of manufacture | Product code | Serial number | Medical device | Note instructions for use eIFU www.meyco.ch | Store in a dry place | Fragile, handle with care | Do not use if packaging is damaged | CE class I _r |
|--------------|---------------------|--------------|---------------|----------------|--|----------------------|---------------------------|------------------------------------|-------------------------|
| | | REF | SN | MD | | | | | CE 1250 |

Description

These diamond knives are reusable surgical instruments comprising a handpiece and a diamond blade. The monocrystalline diamond blades are 100µm to 200µm thick. The handpiece is made of titanium.

Intended use

The product is a surgical cutting instrument (scalpel). During surgery, it comes into contact with human tissue. The handheld product may be used in a variety of procedural steps and techniques in surgery. It is a reusable surgical instrument.

Indications

Diamond knives are used in various surgical procedures, such as

- Ophthalmic surgery
- Refractive surgery (including limbal relaxing incisions).
- Cataract surgery (phacoemulsification, paracentesis, clear cornea incisions)
- Glaucoma surgery (including scleral flap for trabeculectomy)
- Keratoconus (implantation of Kerarings and Intacs)
- Retinal surgery (subretinal membrane transection, vitrectomy)
- Microsurgery

Contraindications

- Cutting of hard tissues, such as bones.

Patient group

No restriction

Users

The diamond knife may only be used by healthcare professionals (ophthalmologists, surgeons). The product is a surgical knife and, as such, requires no product-specific training.

Residual risks / side effects

Risk of cuts if handled improperly.

No other side effects or residual risks are known.

Warnings / precautions

- Read instructions before use.
- The diamond knife is shipped non-sterile. The diamond knife must be thoroughly cleaned and sterilised before use. It must also be cleaned and subsequently sterilized after each

subsequent use. For instructions on cleaning and sterilising the diamond knife, see the sections on **Cleaning and Sterilisation**.

- Never let the diamond blades come into contact with other instruments. Special care must thus be taken during surgery to ensure that the blade does not come into contact with instruments such as forceps or phaco tips.
- **Never store** the diamond knives **uncleaned**.
- Hand the diamond knife to the surgeon with the blade retracted and have him/her return it the same way. **Do not rest** the diamond knife **open** on the instrument tray.
- Before each use, always ensure that the diamond blade is undamaged. Do not use damaged diamond knives.
- **After each use, ensure that the blade of the diamond knife has been retracted into the handpiece with the pressure mechanism or bayonet catch.**

Use / service life / maintenance

Only use the sterilised diamond knife in a correctly prepared sterile setting.

The diamond knife can be reprocessed up to 500 times. Following this, it must be returned to the manufacturer for inspection.

Specific notes on usage

- **Diamond knife with micrometer (ME400 – ME499)**
Adjusting the micrometer barrel:
 1. Rotate the micrometer barrel to the second marker ring
 2. Hold the knife vertically and depress the pressure mechanism
 3. Turn the micrometer barrel to zero – the tip of the diamond will now be flush with the guard rail.
 4. Set the desired cutting depth – 1 revolution = 0.50 mm
We recommend verifying the cutting depth at each use. The cutting depth can be determined with gauge ME409 (for diamond knives ME400, ME401) or ME459 (for diamond knives ME410, ME411, ME420, ME425, ME450, ME451).

Optional accessories for diamond knives with adjustable cutting depth

Gauges

- ME409
- ME459

These gauges can be used in general to verify that the depth setting is working as intended. The gauges cannot be cleaned and sterilised as specified below. Hence, they cannot be employed in sterile conditions.

Instructions for reprocessing

Cleaning

Automated and manual cleaning of the product was performed and validated by SMP GmbH, D-72072, on behalf of Anton Meyer & Co AG. Validation was based on the following criteria:

| No. | Criteria | Accepted threshold value | References |
|-----|---|--------------------------|---|
| 1 | No visible residue on the knife after completion of the cleaning process. | N/A | EN ISO 15883-1:2014 / RKI guideline: 2012 AAMI TIR 30: 2011 /Guideline of DGKH, DGSV, AKI: 2014 |
| 2 | The amount of proteins per sample is below the following value: | <100 µg | RKI guideline: 2012 |
| 3 | The amount of proteins per sample is below the following value: | <200 µg | EN ISO 15883-1:2014 |
| 4 | The amount of proteins/cm ² per sample is below the following value: | < 6.4 µg/cm ² | AAMI TIR 30: 2011 Alfa <i>et al</i> AJIC 1999 |
| 5 | The amount of proteins/cm ² per sample is below the following value: | < 3.0 µg/cm ² | Guideline of DGKH, DGSV, AKI: 2014 |
| 6 | The amount of haemoglobin/cm ² (blood pigment) per sample is below the following value | < 2.2 µg/cm ² | AAMI TIR 30: 2011 Alfa <i>et al</i> AJIC 1999 |
| 7 | The amount of radioactivity per sample is below the following value | < 5 cps | AAMI TIR 11011010605: Validation of the acceptance criteria of the Radionuclide Method as performed by SMP GmbH |

Automated cleaning


Automated cleaning must be carried out using a washer-disinfector (WD) complying with EN ISO/ANSI AAMI ST15883 that is generally acknowledged as effective (e.g., bearing the CE marking according to EN ISO 15883 or approved/supervised/registered by the DGHM or FDA).

1. Pre-cleaning:

- Rinse the instrument under running cold water until there are no more residues.
- Preclean the handpiece under running water with a nylon brush.
- Never brush the diamond blade itself! In addition, never allow it to come into contact with any other object.

2. Loading the cleaning carrier:

- Place the diamond knives with the blade retracted into the metal cleaning carrier and secure them with silicone clips (see figure below).

3. Automated cleaning:

Use the following programme for automated cleaning:

- Preclean for 2 min with cold tap water
- Clean for 5 min with 55°C tap water and 0.5% neodisher MediClean forte cleaning solution (Dr. Weigert, D-Hamburg) or equivalent.
- Empty
- Rinse for 3 min. with cold deionised water
- Empty
- Rinse for 2 min. with cold deionised water
- Empty

4. Disinfection

If possible, select an approved WD programme for thermal disinfection (A0 value ≥ 3000 or – for older devices – at least 5 min at 90°C/194°F; chemical disinfection runs the risk that residual disinfectant may remain on the instruments).

5. Drying

Dry the instruments with an approved WD programme. To prevent re-contamination, ensure only filtered air (free of oil, micro-organisms and other particles) is used.

Manual cleaning

1. Clean the handpiece with a soft nylon brush under cold tap water (blade must be retracted).
Rinse the diamond blade under cold tap water for 10 s or until no visible residue remains. Never brush the diamond blade itself! In addition, never allow it to come into contact with any other object.
2. Immerse the diamond knife for 10 min in 0.5% neodisher MediClean forte cleaning solution (Dr. Weigert, D-Hamburg) at room temperature (20°C ± 2°C).
3. Rinse the knife with deionised water (room temperature 20°C ± 2°C) for at least 1 minute.

Caution:

- **Never allow the diamond blade to come into contact with other instruments, as this may damage the blade.**
- After cleaning, check that the instrument is functioning properly. Then immediately place the diamond knife in the sterilisation box or another appropriate container.
- Automated cleaning is preferable to manual cleaning.

Sterilisation

- The diamond knife must be sterilised before each use.
- The diamond knife must be sterilised with the diamond blade retracted. Thus be sure that the blade is fully within the titanium handpiece.
- The sterilisation trays must be lined with protective materials (knobbed mats).

| Sterilisation method | Sterilisation type | Sample load | Temperature | Sterilisation duration in minutes |
|----------------------|-----------------------|--|---------------|-----------------------------------|
| Steam | Gravitation procedure | Sterile barrier system- Double package: | 134°C (273°F) | 15 minutes |
| Steam | Pre-vacuum procedure | STERKING – R43 type inside. STERKING – R43 type outside. | 134°C (273°F) | 5 minutes |

The sterilisation methods described above were tested and validated according to EN ISO 14937.

Caution:

- The material used to manufacture diamond blades is extremely hard but also highly fragile. Handle the ultra-sharp blades with care to avoid damaging them.
- Avoid physical contact with the blade at all cost.

Storage

There are no specific instructions on storage conditions. Diamond knives should be stored in a dry environment.

Disposal:

Observe country-specific laws and regulations for disposal.

To avoid injury, dispose of diamond knives in a protective disposal container designed for sharp instruments.

Reporting incidents and adverse events

Report all product-related incidents or adverse events immediately to the manufacturer or to the appropriate competent authorities.

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ANTON MEYER & CO. AG