

The limbal relaxing incision knives



ME-435



ME-700, ME-701, ME-704, ME-709




ME-707




ME-700, ME-701, ME-704, ME-709

setting at 0.50 mm: ME-709:
setting at 0.55 mm: ME-704
setting at 0.60 mm: ME-700
setting at 0.65 mm: ME-701




ME-435

micrometer knife with one side footplate for better control



ME-707

settings at 0.50, 0.55 and 0.60 mm



Swiss Diamond Technology
www.meyco.ch



ME-700, ME-701, ME-704, ME-709

preset step knife ME-700

- 1 setting at 0.60 mm
- single sided footplate for better control
- 6 facets 1.0 mm blade
- all titanium handle
- also available with one setting at
0.50 order no.: ME-709
0.55 order no.: ME-704
0.65 order no.: ME-701

Cutt. depth	B	B1	L	D	A	Order No.
0.50 mm	1.00	0.40	2.00	0.20	40°	ME-709
0.55 mm	1.00	0.40	2.00	0.20	40°	ME-704
0.60 mm	1.00	0.40	2.00	0.20	40°	ME-700
0.65 mm	1.00	0.40	2.00	0.20	40°	ME-701



ME-435

micrometer knife ME-435

- micrometer of the highest precision, with an easy legible, black scale
- single sided footplate for better control
- specially designed for LRI
- 6 facets 1.0 mm blade

B	B1	L	D	A	Order No.
1.00	0.40	3.00	0.10	20°	ME-435



ME-707

preset step knife ME-707

- 3 settings at 0.50, 0.55 and 0.60 mm
- single sided footplate for better control
- 6 facets 1.0 mm blade
- all titanium handle

Cutt. depth	B	B1	L	D	A	Order No.
0.50 mm	1.00	0.40	3.00	0.10	20°	ME-707
0.55 mm						
0.60 mm						



Swiss Diamond Technology
www.meyco.ch