



MEYCO SWISS DIAMOND KNIVES

MICS the adequate diamond knives for micro incision coaxial surgery



the mini navigator (3D)




the miniature trapezoid knife



clear corneal knife with parallel sides




- the new meyco diamond knives to proceed micro incision coaxial surgery
- the swiss monocrystalline diamond knives for unlimited time of use for perfect clear corneal and scleral cuts - a contribution to avoid astigmatism




the 3D knives (navigator)

B	L	D	A	order n°
1.30	4.00	0.20	80°	ME-750
1.40	4.00	0.20	80°	ME-751
1.50	4.00	0.20	80°	ME-752
1.60	4.00	0.20	80°	ME-753
1.70	4.00	0.20	80°	ME-754
1.80	4.00	0.20	80°	ME-755
1.90	4.00	0.20	80°	ME-756
2.00	4.00	0.20	80°	ME-757
2.10	4.00	0.20	80°	ME-758
2.20	4.00	0.20	80°	ME-759
2.30	4.00	0.20	80°	ME-760
2.40	4.00	0.20	80°	ME-761



clear corneal knives with trapezoid sides

B	B1	L	D	A	order n°
0.90	1.20	4.0	0.20	90°	ME-770
1.00	1.30	4.0	0.20	90°	ME-771
1.10	1.40	4.0	0.20	90°	ME-772
1.20	1.50	4.0	0.20	90°	ME-776
1.30	1.60	4.0	0.20	90°	ME-773
1.40	1.60	5.0	0.20	44°	ME-267
1.40	1.70	4.0	0.20	90°	ME-777
1.50	1.80	4.0	0.20	90°	ME-774
1.60	1.90	4.0	0.20	90°	ME-778
1.70	2.00	4.0	0.20	90°	ME-779
1.80	2.10	4.0	0.20	90°	ME-775
2.00	2.30	4.0	0.20	90°	ME-781
2.10	2.40	4.0	0.20	90°	ME-782



clear corneal knives with parallel sides

B	L	D	A	order n°
1.30	4.00	0.20	90°	ME-730
1.40	4.00	0.20	90°	ME-731
1.50	4.00	0.20	90°	ME-732
1.60	4.00	0.20	90°	ME-733
1.70	4.00	0.20	90°	ME-734
1.80	4.00	0.20	90°	ME-735
1.90	4.00	0.20	90°	ME-736
2.00	4.00	0.20	90°	ME-737
2.10	4.00	0.20	90°	ME-738
2.20	4.00	0.20	90°	ME-228
2.30	4.00	0.20	90°	ME-229



Swiss Diamond Technology
www.meyco.ch